



**CHOKING MANIKIN CARE
& INSTRUCTION MANUAL**

The Simulaid's Choking Manikins are designed to provide students in first aid classes the opportunity to practice clearing an obstructed airway.

Each manikin is constructed of durable vinyl plastic. A rib cage, sternum and xiphoid process provide the necessary landmarks for practicing the obstructed airway maneuver. Abdominal thrusts and chest thrusts will cause the expulsion of the obstruction.

HOW TO SET UP THE MANIKIN

Select a choking option and drape it over the neck of the manikin. The round ball is the least difficult to expel. The simulated "piece of meat" increases the volume of the obstruction and the number of thrusts needed to free the airway. Use one option at a time so the cords do not become entangled. Note: Infant Choker size will only include and utilize the round ball.

CARE AND MAINTENANCE

Since the manikin should only be used to practice obstructed airway maneuvers, it will only be necessary to clean the exterior surfaces. The head and torso should be washed with warm, soapy water, rinsed thoroughly with clear water and dried.

Caution: Do not allow water to enter the internal parts of the foam filled manikin.

The carry bag included with your product is first quality, colorfast nylon. The physical properties of your new manikin will, occasionally, pick up a trace of color from your bag. This will NOT affect the use and viability of your manikin or your carry bag.

PLEASE STORE MANIKIN IN UPRIGHT POSITION.

No. 1602 Adult Size Torso with Carry Bag

Size: 31" x 18" x 10" Sh. wt. 16 lbs.

No. 1630 Obese Mature Torso with Carry Bag

Size: 30" x 21" x 13" Sh. wt. 20 lbs

No. 1615 Adolescent Size Torso with Carry Bag

Size: 24" x 14" x 9" Sh. wt. 9 lbs.

No. 1620 Child Size Torso with Carry Bag

Size: 17" x 10" x 7" Sh. wt. 5 lbs.

No. 1640 Infant Size Full Body with Carry Bag

RETURN POLICY

Should it be necessary to return an item, contact our Customer Service Department to obtain an RGA Number. Please refer to your invoice number when phoning in your request for returning merchandise. Should you have any questions or need further information on any product we manufacture, call or write our Customer Service Department.

WARRANTY

Simulaids warrants this product to be free from any defect in materials and/or workmanship for a period of three years from the date of purchase, as evidenced by the date of invoice when the product was shipped to the end user. This warranty expressly does not cover abuse, accidental or purposeful damage, or any form of modification to the product. Simulaids reserves the right to either repair or replace affected parts or the entire unit, at their sole discretion, after investigating and reviewing the actual product and the damage. In most instances, a digital photo of the product in question showing the damage will help qualify a product for return to the factory. At no time will any product be accepted at the plant without proper return authorization issued by Simulaids.

Freight and Shipping charges are the sole responsibility of the end user. No product will be received with shipping charges due. Any product considered for warranty work must be identified by serial number and invoice number from the agency through whom the product was purchased. Without this information the product will not receive a return authorization number as required above.



PO Box 1289 - 16 Simulaids Drive

Saugerties New York 12477

800-431-4310 or 845-679-2475

Fax 845-679-8996

e-mail info@simulaids.com

www.simulaids.com



 Simulaids



www.simulaids.com