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[EN] PRECAUTIONS

Please read the instructions carefully before using this device. To reduce the risk of accidents, take the following precautions:

1. Warning: The IEC 60417 - 5582 symbol indicates that this device should not be used near bathtubs, showers, sinks or other containers containing water.
2. Do not place or drop it in water or other liquids.
3. If the device falls into water, unplug it immediately. Do not put your hand in the water.
4. Keep the device away from water.
5. The device must never be left unattended while plugged in.
6. Use this device only for the purpose for which it was designed, as indicated in this manual. Do not use accessories not recommended by the manufacturer.
7. The appliance must never be left unattended while plugged in.
8. This appliance can be used by children aged 12 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been supervised or instructed on how to use the appliance safely and understand the dangers involved.

9. Never operate this appliance if the cord or plug is damaged, if it is not working properly, if it has been dropped or if it is in any way damaged or has fallen into water. Take the appliance to the nearest authorised service centre for examination or repair.

10. Keep the cable away from heat or hot surfaces. Do not wrap the cable around the appliance.

11. Do not use an extension cable with this appliance.

12. Store the appliance in a dry, dust-free place away from direct sunlight.

13. Do not attempt to open the control compartment under any circumstances and never attempt to repair the equipment yourself.

14. Do not look directly at the LEDs when they are lit.

15. Keep cables away from heat or hot surfaces.

Do not wrap cables around the device.

16. Never use if you are drowsy or half asleep.

17. Store the product in a dry environment at room temperature, avoiding dust, direct sunlight or humid environments.

18. This product contains a built-in rechargeable battery that cannot be replaced or removed. Do not throw this device into a fire, charge it, use it, or place it in high temperatures.

19. Avoid looking directly at the light source of the mask or looking at the light for a long time. Keep your eyes closed during use. If you feel that the light is too strong, cover your eyes with cotton wool or an eye mask.

20. Because the mask material is plastic, the opened packaging may have an unpleasant odour, which is normal, but the odour will quickly disappear.

21. Clean the mask after use with a soft, slightly damp cloth. The elastic band can be washed with water.

22. Do not rub the mask with corrosive chemical solvents.

23. Avoid exposure to sunlight.

24. Consult your doctor before using this device. Do not use this device with:

- Pregnant women
- Children
- People with epilepsy

- Patients with thyroid problems or who are taking thyroid medication should not use it
- People with light allergies or who are taking medication that can cause light allergies
- People with light-sensitive skin.
- People with skin wounds or any adverse condition.
- People with eye diseases.

25. This device is equipped with a rechargeable battery that cannot be replaced by the user. In case of malfunction, contact technical support.

26. Do not disassemble, modify or remove the battery during use or the device's service life. Do not insert foreign objects or immerse the battery in water or other liquids. These actions may cause leakage, overheating, fire or explosion.

27. When the device reaches the end of its useful life, the battery must be removed from the device before disposal, for proper management as hazardous waste. Ensure that the device is completely turned off.

28. The battery must be disposed of safely in accordance with local regulations on electronic waste.

29. If you are unable or unwilling to remove the battery yourself to avoid risks, take the equipment to an authorised recycling centre, where qualified personnel can safely remove and dispose of the battery.

30. Ensure that the battery is fully charged at least once every three months, otherwise the battery may be damaged.

31. Keep these instructions for future reference.

[EN] MODO DE USO

LIGHT SENSITIVITY TEST:

IMPORTANT: Before using this product for the first time, perform a light sensitivity test:

Place the device on your arm, turn it on and select a treatment time of 10 minutes. If your skin is red or itchy after this time, it may indicate that you have photosensitive skin. We recommend that you stop using this device until you have consulted a dermatologist before using it again.

INSTRUCTIONS FOR USE:

1. Before using the mask, the battery must be recharged for it to work. The LED mask will not work while it is plugged in.

2. Attach the ends of the elastic Velcro straps to the LED mask using the holes on the left and right sides of the mask.

3. Cleanse your face and allow the water to dry naturally or use an absorbent towel to gently dry your face (do not rub your face vigorously, as this will cause your pores to become clogged or your skin to become rough and dehydrated).

4. To turn on the device, press the touch button on the front with your finger for 3 seconds. The colour will change once the machine is turned on. The colour will change automatically after the eighth touch.

5. Cover your face with the light mask and find a comfortable and suitable position. Pass the elastic band behind your ears and adjust the tension using the sliders at the ends.

6. We recommend performing the photon beauty treatment during your rest time. Find a comfortable sofa or bed, lie down comfortably, close your eyes, and enjoy the beauty process provided by the photons. Do not read or watch television for long periods of time while using the device.

7. You can use cosmetic products on your skin together with the LED mask to amplify their properties and absorption. Apply them evenly to your face after cleansing. We recommend using red light for best absorption. Generally, we recommend using it for 45 to 30 minutes.

8. Remove the mask, turn it off, and clean it so you can use it again next time.

USE OF COLOURS:

Red: 630NM. Helps skin look firmer and more luminous. Contributes to a more elastic feeling on the face, promotes softness and supports the improvement of the appearance of lines and wrinkles.

Blue: 470NM. Helps to calm imperfections and supports the visual cleansing of the complexion. Promotes the reduction of redness associated with breakouts and helps the skin look more even.

Green: 520NM. Helps balance the appearance of the skin, promotes a feeling of calm and helps regulate excess shine. It also promotes the reduction of facial swelling.

Yellow: 585NM. Promotes a more even and luminous appearance, helps smooth texture and contributes to improving areas with redness. It also promotes a more revitalised appearance.

Purple: 426NM. Combines the benefits of red and blue light. Helps improve the appearance of marks and imperfections while promoting a more balanced look.

Turquoise: 490NM. Helps skin feel fresher and revitalised. Contributes to promoting the natural dynamism of the complexion.

White: Combination of all colours. Promotes deeper illumination of the face, helps improve the appearance of dull skin and contributes to a feeling of firmness and softness.

[EN] ENVIRONMENT

Disposal of old Electrical & Electronic Equipment (Applicable throughout the European Union and other European countries with separate collection programs. This symbol, found on your product or on its packaging, indicates that this product should not be treated as household waste when you wish to dispose of it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent negative consequences to the environment and human health, which could otherwise be caused by inappropriate disposal of this product. The recycling of materials will help to conserve natural resources. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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